(EN) MANUAL - PULSE OXIMETER A3 - INSTRUCTIONS TO USER

IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!

Read the instruction manual carefully before using this device, especially the safety instructions and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

Failure to follow these instructions can cause damage to the device and personal injury. The manufacturer is not responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and damage to the device due to user's negligence of the operation instructions. Intended use: spot checking of oxygen saturation, SpO₂ and pulse rate of an adult. Not for continuous monitoring.

WARNING:

- An uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.
- The device cannot be clipped on the edema and tender tissue.
- The infrared light emitted from the device is harmful to the eyes. The user cannot stare into the light.
- Users fingernail cannot be too long and must be clean.
- Fingernail polish or acrylic nails may result in inaccurate readings.

1. SAFETY

1.1 INSTRUCTIONS FOR SAFE OPERATIONS

- Check the device and all accessories periodically to make sure that there is no visible damage that may affect user's safety and monitoring performance about cables and transducers. It is recommended that the device should be inspected at least once a week. When there is obvious damage, stop using the device.
- Necessary maintenance must be performed by qualified service engineers only. Users are not permitted to maintain it by themselves.
- The device cannot be used together with devices not specified in the user manual. Only the accessory supplied with this device can be used.
- The equipment is calibrated to display functional oxygen saturation.

1.2 WARNINGS

- Explosive hazard - DO NOT use the device in environments with inflammable gas such as some ignitable anesthetic agents.
- DO NOT use the device in an MRI or CT environment.
- DO NOT use this device when the user is allergic to rubber.
- The disposal of scrap instrument, accessories and packings (including battery, plastic bags, foams and paper boxes) should follow the local laws and regulations.
- Please check the packing before use to make sure the device and accessories are totally in accordance with the packing list, or else the device may have the possibility of working abnormally.
- Please do not measure this device with function test paper for the device's related information.

1.3 ATTENTIONS

- Keep the device away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- If the device gets wet, please stop operating.
- When it is carried from cold environment to warm or humid environment, please do not use it immediately.
- DO NOT operate keys on front panel with sharp materials.
- High temperature or high-pressure steam disinfection of the device is not permitted. Refer to the user manual in the relative chapter for instructions of cleaning and disinfection.

- Do not immerge the device into liquid. When it needs cleaning, use a soft cloth, moistened with isopropyl alcohol. Do not spray any liquid on the device.
- When the fingers are too thin or too cold, it will affect the normal measure of the users SpO₂ and pulse rate, please clip a finger such as the thumb or middle finger deeply enough into the probe.
- Do not use the device on infants or neonatal patients.
- The device is suitable for children above 14 years old and adults.
- The device may not work for all users. If you are unable to achieve stable readings, discontinue use.
- The update period of data is less than 5 seconds, which is changeable according to different individual pulse rate.
- The waveform is normalized. Please read the measured value when the waveform on screen is equally and steady going, here this measured value is optimal value and the waveform at the moment is the standard one.
- If some abnormal conditions appear on the screen during test process, pull out the finger and reinsert to restore normal use.
- The device has normal useful life for three years since the first electrified use.
- Batteries must be removed if the device will be stored for more than one month, or else batteries may leak.
- A flexible circuit connects the two parts of the device. Do not twist or pull on the connection.

2. OVERVIEW

2.1 CLASSIFICATION

Class II B, (MDD93/42/EEC IX Rule 10)

2.2 FEATURES

- The total weight is +/- 50 g including batteries.
- Power consumption of the product is low. With two AAA batteries the device can operate continuously for 50 hours.
- The device will enter standby mode when no signal is in the product within 8 seconds.
- Display direction can be changed, easy to view.

2.3 MAJOR APPLICATIONS AND SCOPE OF APPLICATION

The device can be used to measure human hemoglobin saturation and pulse rate through the finger and indicate the pulse intensity by the bar-display.

The device is not suitable for use in continuous supervision for patients.

The problem of overrating would emerge when the patient is suffering from toxicosis which caused by carbon monoxide, the device is not recommended to be used under this circumstance.

2.4 ENVIRONMENT REQUIREMENTS

Storage environment

- a. Temperature: -20° C ~ +55° C
- b. Relative humidity: 10% ~ 93%
- c. Atmospheric pressure: 22 kPa ~ 107.4 kPa

Operating environment


- a. Temperature: 0° C ~ 40° C
- b. Relative humidity: 15% ~ 80%
- c. Atmospheric pressure: 59 kPa ~ 107.4 kPa

3. PRINCIPLE AND CAUTION

3.1 PRINCIPLE OF MEASUREMENT

Principle of the device is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin (Hb) and oxyhemoglobin (HbO₂) in glow & near-infrared zones. Operation principle of the instrument is: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning & Recording Technology, so that two beams of different wavelength of lights can be focused onto human nail tip through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.

Figure 1. Operating principle

3.2 CAUTION

1. The finger should be placed properly (see the attached illustration of this manual, Figure 5), or else it may cause inaccurate measurements.
2. The SpO₂ sensor and photoelectric receiving tube should be arranged in a way with the subject's arteriole in a position there between.
3. The SpO₂ sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
4. Make sure the optical path is free from any optical obstacles like rubberized fabric.
5. Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight etc.
6. Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.
7. **Pulse rate and SpO₂ alarm**

When a certain physiological parameter exceeds the set alarm range, a physiological alarm is triggered, and the parameter on the screen will change color, reminding the user to pay attention to the measured parameter.

SpO₂ parameter values:

Color	Range
Green	85-99
Orange	Parameter value is outside the 'green range'

Pulse Rate value:

Color	Range
Green	50-120
Orange	Parameter value is outside the 'green range'

3.3 CLINICAL RESTRICTIONS

1. As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.

2. For those with a substantial amount of staining dilution drug (such as methylene blue, indigo green and acid indigo blue), or carbon monoxide hemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic hemoglobin, and some with icterus problem, the SpO₂ determination by this monitor may be inaccurate.
3. Drugs like dopamine, procaine, prilocaine, lidocaine and butacaine may also be a major factor blamed for serious error of SpO₂ measure.
4. As the SpO₂ value serves as a reference value for judgement of anemic anoxia and toxic anoxia, some patients with serious anemia may also report good SpO₂ measurement.

4. TECHNICAL SPECIFICATIONS

1. **Display format:** TFT Display; **SpO₂ measuring range:** 35% ~ 99%; **Pulse rate measuring range:** 30 bpm ~ 240 bpm;
2. **Power requirements:** 2x 1.5 V AAA battery
3. **Power consumption:** 150 mW
4. **Resolution:** 1% for SpO₂ and 1 bpm for pulse rate; perfusion index: 0.1%.
5. **Measurement accuracy:** ±2% in stage of 70% ~ 100% SpO₂ and meaningless when stage being smaller than 70%. ±2 bpm during the pulse rate range of 30 bpm ~ 99 bpm and ±2% during the pulse rate range of 100 bpm ~ 250 bpm.
6. **Resistance to surrounding light:** The deviation between the value measured in the condition of man-made light or indoor natural light and that of darkroom is less than ±1%.
7. **It is equipped with a function switch:** The device can enter standby mode in case of no finger is inserted in the Oximeter within 8 seconds.
8. **Optical sensor:** Red light (wavelength is 660 nm, 6.65 mW) Infrared (wavelength is 880 nm, 6.75 mW)

5. INSTALLATION

5.1 VIEW OF THE FRONT PANEL

Figure 2. Front view



Figure 3. Batteries installation

5.2 BATTERY

- Step 1. Insert 2x AAA size batteries properly in the right direction.
- Step 2. Replace the cover.

Please take care when you insert the batteries. The improper insertion may damage the device.

6. OPERATING GUIDE

Figure 5. Put the finger in position

1. Clip the finger in between the rubber cushions (make sure the finger is in the right position).
2. Press the switch button on front panel to exit standby mode.
3. Do not move the finger and keep the user at ease during the process.
4. Read the information directly from the screen display.
5. The device could change display direction according to the handing direction.

Fingernails and the luminescent tube should be on the same side. Measured result. The values measured with this device are not suitable for making diagnosis. The oxygen saturation (SpO_2) is a term referring to the concentration of oxygen in blood. The value normally lies between 90% and 96%.

7. REPAIRING AND MAINTENANCE

- Change the batteries when the low voltage is displayed on the screen.
- Clean the surface of the device before using. With a soft cloth, moistened with isopropyl alcohol, let it dry in air or clean it by dry clean fabric. To prevent a cross infection from the previous use.
- Please take out the batteries if the device is not in use for a long time.
- The best storage environment of the device is -20°C to $+55^\circ\text{C}$ ambient temperature and not higher than 93% relative humidity.

High-pressure sterilization cannot be used on the device. Do not immerse the device in liquid. It is recommended that the device should be kept in a dry environment. Humidity may reduce the useful life of the device, or even damage it.

8. TROUBLESHOOTING

Trouble	Possible Reason	Solution
The SpO_2 and pulse rate cannot be displayed normally	1. The finger is not positioned properly. 2. The users SpO_2 is too low to be detected.	1. Place the finger properly and try again. 2. Try again.
The SpO_2 and pulse rate are not displayed stably	1. The finger is not placed inside deep enough. 2. The finger is shaking or the user is moving.	1. Place the finger properly and try again. 2. Let the patient keep calm.
The device cannot be turned on	1. The batteries are drained or almost drained. 2. The batteries are not inserted properly. 3. The malfunction of the device.	1. Change batteries. 2. Reinstall batteries. 3. Please contact the local service center.
The display is off suddenly	1. The device will enter standby mode when no signal is in the product within 5 seconds. 2. The batteries are almost drained.	1. Normal. 2. Change batteries.

9. KEY OF SYMBOLS

Symbol	Description
	Type BF
	Refer to instruction manual/booklet
$\%\text{SpO}_2$	The pulse oxygen saturation(%)
PRbpm	Pulse rate (bpm)
PI	Perfusion Index (%)
	The battery voltage indication is deficient (change the battery in time avoiding the inexact measure)
	1. No finger inserted 2. An indicator of signal inadequacy
	Battery positive electrode
	Battery cathode
	1. Change brightness of the screen. 2. Exit standby mode.
SN	Serial number
	WEEE (2002/96/EC)
IP22	International protection
	This item is in according to the Medical Device Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community.
	Manufacturer
	Date of manufacture
	Storage and Transport Temperature limitation
	Storage and Transport Humidity limitation

10. FUNCTION SPECIFICATION

Display information	Display mode
Pulse Oxygen Saturation (SpO ₂)	TFT
Pulse rate (PR)	TFT
SpO ₂ parameter specification	
Measuring range	0% ~ 100%, (the resolution is 1%)
Accuracy	70% ~ 100%: ±2% Below 70% unspecified
Optical sensor	Red light (wavelength is 660 nm) Infrared (wavelength is 880 nm)
Pulse parameter specification	
Measuring range	30 bpm ~ 250 bpm (the resolution is 1 bpm)
Accuracy	±2 bpm or ±2% select larger
Pulse intensity	
Range	Continuous bar-graph display, the higher display indicates the stronger pulse.
Battery requirement	
1.5 V (2x AAA size) batteries	
Battery useful life	
Two batteries can work continually for 20 hours	
Dimensions and weight	
Dimensions	63 × 41 × 31 mm

Guidance and manufacturer's declaration - electromagnetic emission		
The pulse oximeter is tended for use in the electromagnetic environment specified below. The customer of the user of the pulse oximeter should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The pulse oximeter uses RF energy only for their internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The pulse oximeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emission IEC 61000-3-3	Not applicable	

Shenzhen Hexin Zondan Medical Equipment Co.,Ltd
Floor 14, Block D,Dianlian Technology Building, the Crossing
between South Circle Road and South Fuli Road, Guangming
District, Shenzhen, CHINA
Post Code:518067
Tel: +86-0755-26865970, 26850278

Shanghai International Holding Corp. GmbH(Europe)
Eiffestrasse 80, 20537 Hamburg Germany

Mascot Online BV
Veluwezoom 32
1327 AH Almere
The Netherlands
www.mascot-online.nl

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