WellHealth

OWNER'S MANUAL

Wellhealth mini massage gun

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Bullet head/ For all body parts

Instructions for use

1. Power On/Off button (Long press to power on/Long press to power off)
2. Press the power button to change the mode
3. To remove the detachable massage heads, please just pull it out with force.
4. Please use a 5V charging device to charge it.
5. Please charge it for three hours before the first use. Please connect it to a 5V adaptor with a cable for charging.

Safety instructions

To reduce the risk of electric shock, fire, and personal injury, you must follow instructions when using this product.

For adults only

- This product can only be used on dry and clean skin, don't use it on clothes. Gently move it on your skin. Keep massaging for about 60 seconds at each part of your body.
- Only using on soft tissues of your body shall not cause pain and discomfort. Please do not use it on head or any hard or bony part of your body
- Only use the massage head generating the best effect.
- Bruises may occur under moderate or normal pressure. Please always check the area for physical therapy.
- Please discontinue the massage at the first sign of pain or discomfort.
- Please always keep your fingers, hair or other body parts away from the axle and back of the massage head, as they may be squeezed.
- Don't place any object in the vent of the instrument.

- Do not soak it in water or let water get into the vent hole of the instrument.
Don't drop the product or use it improperly. Only 5V DC charger can be used to charge the instrument.
- Please check the instrument and its batteries carefully before each use. Do not change or modify the instrument in any way.
- Never leave it being operated or charged unwatched.

Indications and contraindications

Deep tissue massage of the massage gun can provide pain relief for joint sprains, bursitis, muscle cramps, neuralgia, and other musculoskeletal disorders in a short time. Almost all soft tissue disorders are associated with trauma, overuse of muscles, or bad posture. This massage gun can be used to treat most scar tissue associated with these disorders. Please avoid excessive pressure and prolonged physical therapy during use. Please also avoid sensitive parts such as the head, face and cervical spine. Stop the physical therapy in case of inflammation, swelling or increased pain.

Ongoing research on the effects of vibration training on specific medical conditions is likely to reduce the types of contraindications described above, and practical experience indicates that in many cases, the integration of vibration training with physical therapy regimens is also a care treatment recommendation, which must be completed in the presence of a doctor, specialist or physical therapist.

Indications	Contraindications
<ul style="list-style-type: none"> » Pain and cramps from muscle injuries, sprains, strains » Relax the thickened connective tissue and fascia » Reduce the buildup of lactic acid in muscles » Increase the range of joint motion » Relieve muscle fatigue 	<ul style="list-style-type: none"> » Aneurysms, bleeding, use of blood thinner » Heart disease with a pacemaker or defibrillator » Pregnancy, cancer » Parts close to bone and with less muscle coverage, such as tibia, instep, and the back of the hand

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