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



Product Instruction

Device Introduction

Device Name	IPL Ice-cool Hair Remover	Device Size	136.5x77x47.5mm
Power Adapter	100~240V	Illumination Area	3cm ²
Rated Input	DC12V $\approx$ 3A	Lifespan	999,900 flashing
Max Output	36W	Net Weight	Around 249g

Working Principle

- ▶ Under the illumination of impulse light with specific wave band and pulse width, the melanin in the skin absorbs energy and transfers the transient heat to the follicle. When the hair follicles reach certain temperature, they will enter the dormancy period, which makes the hair stop growing by stages and fall off in the end. The ideal effect can be achieved by using the device for around 8 weeks or more, decided according to the physique of different people. The built-in filter will filtrate out ultraviolet light and other useless light, which makes the device ideal for skin care.
- ▶ Ice-cool function can relieve the pain of burning sensation from higher gear. It takes around 10 seconds for the cooling system to cool up.
(Note: even for the same power level, burning sensation may vary from person to person)

Skin Color Chart

 White or alike ✓

 Beige ✓

 Light brown ✓

 Medium ✓

 Dark brown ✓


 brown ✗

Pre-use Notes


- ▶ Do not charge the device since it's plug in. Please unplug the power cord after operation and well store them in the shade.
- ▶ Please shave the hair off first before operation and please go with shaving cream. Otherwise, the skin would be damaged under dry condition.
- ▶ Due to different kinds of skin, for first time user, please run a trial test on small area then use it on a large scale if without adverse reactions in 12 hours.

- ▶ Do not try high gear directly for the first time. Go with the lowest one then adjust the gear accordingly.
- ▶ During operation, do not make illustration on the same area repeatedly. Try to irradiate every inch of skin evenly.
- ▶ Ice-cool function can efficiently relieve the pain (if any, feeling may vary), thus improving your experience of our product.
- ▶ To achieve the best effect, please keep using this device for around 8 weeks or more. At that point, hair would be thinner and start to fall off. After that, remember to reuse it once within 15-30 days to maintain the effect, depending on how fast the hair grows back.
- ▶ Please shave the hair off first if it's longer than 1mm. Otherwise, the long hair may absorb the heat and get burned, which would produce an odor. Sometimes it may also cause skin redness and swelling. This will greatly affect your experience of using our product.
- ▶ Keep distance from children during operation since the flash may hurt the eyes without protective goggles.

Operating Guide

1. Shave the hair off first then make sure the skin is dry and clean.


2. Plug in the power cord, then long press the power button (on the side) to turn on the device, switch the gear by short pressing the power button; short press the ice-cool button to switch on and off the ice-cool function (Flashing and Ice-cool function both on); long press the ice-cool button for 3 seconds to enter "Ice-cool Mode" (Flashing off, only Ice-cool function on).

3. Long press the flashing button to enter "Auto Mode"; long press again to exit this mode.


4. Protective goggles must be put on to protect the eyes when operating.

5. During operation, short press the flashing button to start. No need to press button if in "Auto Mode". The device only works when the flashing outlet is fully covered by skin 90°. At that point, the blue light between the power button and ice-cool button will keep blinking. If not, please adjust the position of flashing outlet until the blue light blinks.

Suggestion for duration

Notice: Specific duration may vary from person to person. The above data are for reference only.

Hair Removal Cycle

Week 1-2	Three times a week	Hair grows slower	
Week 3-4	Twice a week	Hair becomes thinner and softer. Some start to fall off	
Week 4-8 or above	Once within 15 to 30 days	The pores start to shrink and hair stops growing	

Attention: the data above may vary from person to person!

Not Applicable

- ▶ Female under pregnancy, period and lactation, light-sensitive person, juvenile and people with cicatrix are not allowed to use.
- ▶ People with heart disease, diabetes, hypertension, malignant tumor, folliculitis, infectious skin disease, burned skin, varicosity, anapetia, herpes simplex, skin lesions, hematocus and anticoagulation disorder are prohibited to use this device.
- ▶ If other methods of hair removal are applied lately, please do not use this device until hair grows back.
- ▶ If products with whitening element or Vitamin A are applied, please wait one week from the day of the application.
- ▶ This device can't be used with essential oil products.
- ▶ Patients taking aspirin, immunosuppressive drugs, oversensitive drugs (such as tetracycline, sulfonamides, etc.) are not recommended to use this device for the time being.
- ▶ People with insolation or sunburn within last two weeks must wait until the skin recovers.
- ▶ It is forbidden to apply the device on sensitive parts or private parts (such as eye, lip, nipple, vagina, etc.).
- ▶ People with polytrichia or other serious complications are not recommended.
- ▶ People with sensitive skin should pay attention to the ice-cool function. Stop using this product if you don't feel comfortable and consult your doctor for advices.

Attention

- ▶ In winter, people with dry skin should pay attention to daily skin care like maintaining moisture.
- ▶ Please wear goggles when operating. Do not directly stare at the flashing light, which would do great damages to the eyes.
- ▶ Though slight redness and itching of the skin rarely occurs, it is normal phenomenon. If so, wait for the skin to recover before reusing the device.

- Do not spill alcohol directly on the flashing outlet when cleansing. Please use cotton swabs or dust-free cloth dipped in a small amount of medical alcohol (or use alcohol cotton) to clean the stain.
- Avoid using gel cosmetics or liquids when operating the devices so as not to damage your skin or the flashing outlet.
- Avoid eating food like spinach, rape or mustard before operation.
- Adjust the gear accordingly, if the skin allows, the higher the gear, the better the effect. People with darker skin color should be careful when choosing high gear; otherwise it is easily to get the skin burned.
- After operation, do not touch hot water within 6 hours and pay attention to skin humidification and sunscreen.
- Avoid ultraviolet and computer radiation after operation and wear sunscreen when going out.
- During operation, the edge of flashing outlet may get stained or black spot because of the residual hair or dirt left on the edge. Although this doesn't affect functions of the device, please keep the flashing outlet clean at all time.
- This device is not water proof, so please do not use it in moisture environment like shower room.
- If any malfunction takes place, send the device back to manufacturer or dealer to confirm whether it can be repaired. Safety of people always comes first under any circumstances.
- Ensure that the device is far from sites with strong magnetic field, moisture, heavy dust, or high temperature.
- Do not make strong impact or knock on this device. It is prohibited to use it with purposes other than designated care.
- Besides the authorized personnel, others shall not open up the device; otherwise there will be a risk of electric shock. It is also strictly forbidden for non-professionals to carry out maintenance without authorization.

Common problems and solutions

Problems	Reasons	Solutions
Unable to turn on the device	The power supply is not connected	Make sure the adapter or power cord is plugged in
The adapter is hot in use	Normal phenomenon	No actions required
Stains found on the edge of flashing outlet	Hair or dirt left on the edge and get burned during operation	Use alcohol cotton to clean up the edge before operation
Effect is not satisfactory	Gear is too low; Operation is not frequent enough	Adjust to suitable gear; Follow the guidance strictly; Make sure the skin is evenly flashed
No illumination	Flashing outlet is not fully covered by skin	Make sure the flashing outlet forms an angle of 90 degrees on the skin
Automatic overheating protection (Error 3)	Self-protection program for the device being overheated	Put the device in a ventilation place to cool down
Illustration causes allergy	Gear is too high; Flashing outlet is damaged	Adjust to suitable gear; Contact customer service for help
Odor comes with illustration	The flashing outlet is covered by dust; There is hair left on skin	Use alcohol cotton to clean up the edge; Make sure the hair is shaved clean; Adjust to the suitable gear

Warranty

The warranty period for this device, including the power adapter and accessories, is one year from the date of receipt. If a quality issue occurs during the warranty period, we will repair or replace the product free of charge. However, damage caused by improper use or incorrect operation by the customer is not covered under the warranty, and repair services may be subject to additional charges accordingly. Please keep the invoice or receipt and warranty card properly for future warranty service reference.

We may not provide a warranty or return service in the following cases:

- ▶ Damages caused by uncontrollable factors such as fire, earthquake, and flood.
- ▶ Damages caused by other force majeure events.
- ▶ Device malfunction or damages caused by customer's fault, misuse, error, or other abnormal conditions.
- ▶ Damages resulting from customer's self-disassembly or improper storage/handling (such as rodent damage, liquid ingress, etc.).
- ▶ Damages resulting from use under unsuitable conditions, including but not limited to exceeding workload limits.
- ▶ We shall not be liable for any indirect or consequential damages arising from non-product-quality issues, including but not limited to loss of business interests, interruption of business operations, or other losses caused by the operation of the device. Liability for accidents shall be determined in accordance with applicable laws.

After-sales service card

In order to provide customer with better after-sales service, please fill out the information below (all required):

Return Purchase ☐ Exchange Purchase ☐ Repair ☐

Warranty Period	One year from receipt	Date of receipt:
Customer Info	Address	
	Name	
	Telephone	
	Order number	
Reasons		

Certificate
Inspector: 

GO DIGITAL, STAY UPDATED

Scan the QR code below to visit our Deki Support Hub for Full User Manual in multiple languages.

By choosing digital resources, you are working together with us to reduce paper waste and protect our planet.

SCAN ME


Offline Support: If you cannot access the digital versions or require a printed manual in a specific language, please contact us at support@dekilife.com. We are committed to ensuring you have all the information you need.

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Upute za rad

1.Shave the hair off first then make sure the skin is dry and clean.

2.Plug in the power cord, then long press the power button (on the side) to turn gear by short pressing the


the ice-cool button to switch on and off the ice-cool function(Flashing and ice-cool function both on); button for 3 seconds to

3.Long press the flashing button to enter "Auto Mode"; long press again to exit this

ñ.Zaštitne naočale moraju be put on to protect the eyes when operating.


5.During operation,short press the flashing button to start. No need to press button if in "Auto Mode".

izlazni signal je Culls

point, the blue light between the power button and ice-cool button will keep blinking. If not, please

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